

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

DNASE2B Fusion Protein



Basic Information

Catalog Number:

Ag18595

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MMARLLRTSFALLFLGLFVGLGAATISCRNEEGKAVD
WFTFYKLPKKQNKESGETGLELYLDSTTRSWRKSEQL
MNDTKSVLGRTLQQLYEAYASKSNNTAYLIYNDGVPK
PVNYSRKYGHTKGLLLWNRVQGFWLIHSIQFPPIPEE
GYDYPPTGRRNGQSGICITFKYNQYEAIDSQLLVCNP
NVYSCSIPATFHQELIHMPQLCTRASSEIPGRLLTTL
QSAQGQKFLHFAKSDSFLDDIFAAWMAQRLKTHLLT
ETWQRKRQELPSNCSLPYHVYNIKAIKLSRHSYFSSYQ
DHAKWCISQKGTKNRWTCIGDLNRSHPHQAFRSGGFI
CTQNWQIYQAFQGLVLYYESCK
(1-361 aa encoded by BC126136)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.