

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RASD1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag18818

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

NKGD R D F Y R E V D Q R E I E Q L V G D D P Q R C A Y F E I S A K K N  
S S L D Q M F R A L F A M A K L P S E M S P D L H R K V S V Q Y C D V L  
H K K A L R N K K L L R A G S G G G G D P G D A F G I V A P F A  
(145-250 aa encoded by BC018041)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.