

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PAX2 Fusion Protein



## Basic Information

**Catalog Number:**

Ag19011

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MDMHCKADPFSAMHPGHGGVNLGGVFVNGRPLP  
DVVRQRIVELAHQGVRPCDISRQLRVSHGCVSKILGR  
YYETGSIKPGVIGGSKPKVATPKVVDKIAEYKRQNPTM  
FAWEIRDLLAEGICDNDTVPSVSSINRIIRTKVQQPF  
HPTPDGAGTGVTAPGHTIVPSTASPPVSSASNDPVG  
SYSINGILGIPRSNGEKRRKREVEVYTDPAHIRGGGGL  
HLVWTLRDVSEGSVPNGDSQSGVDSLKHLRADTFT  
QQQLEALDRVFERPSYDPVFQASEHIKSEQGNEYSLP  
ALTPGLDEVKSSLSASTNPELGSNVSGTQTPVVTGR  
DMASTTLPGYPPHVPTGQGSYPTSLAGMVPGEFS  
GNPYSHQYTAAYNEAWRFSNPALLMPPPGPLPLLP  
MTATSYRGDHIKIQADSFGLHIVPV  
(1-431aa encoded by BC148710)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.