

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PGP Fusion Protein



## Basic Information

**Catalog Number:**

Ag19272

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MAAAEAGGDDARCVRLSAERAQALLADVDTLLFD  
CD  
GVLWRGETAVPGAPEALRALRARGKRLGFITNNS  
SKT  
RAAYAEKLRRLGFGGPAGPGASLEVFGTAYCTAL  
YLR  
QRLAGAPAPKAYVVGSPALAAELEAVGVASVGV  
GPE  
PLQGE  
GPGDWLHAPLEPDVRAVVVGFDPHFSYMKLT  
KALRYLQPGCLLVGTNMDNRLPLENRFIAGTGC  
LV  
RAVEMAAQRQADIIGKPSRFIFDCVSQEYGINP  
ERTV  
MVGDRLDTDILLGATCGLKTLTLTGVSTLGDV  
KNNQE  
SDCVS  
KKKMPDFYVDSIADLLPALQG  
(1-321aa encoded by BC157035)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.