

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

LAMC3-Specific Fusion Protein



Basic Information

Catalog Number:

Ag19464

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

SLSEGGCRPCTCNPAGSLDTC DPRSGRPC KENVEG
NLCDRCRPGTFNLQPHNPAGC S S C F C Y G H S K V C A S T
A Q F Q V H H I L S D F H Q G A E G W W A R S V G G S E H P P Q W S
P N G V L L S P E D E E E L T A P E K F L G D Q R F S Y G Q P L I L T F R V P
P G D S P L P V Q L R L E G T G L A L S L R H S S L S G P Q D A G H P R E V
E L R F H L Q E T S E D V A P P L P P F H F Q R L L A N L T S L R L R V S P G
P S P A G P V F L T E V R L T S A R P G L S P P A S W V E I C S C P T G Y T
G Q F C E S C A P G Y K R E M P Q G G P Y A S C V P C T C N Q H G T C
D P N T G I C V C S H H T E G P S C E R C L P G F Y G N P F A G Q A D D
C Q

(421-753 aa encoded by NM_006059)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.