

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# SNAP47 Fusion Protein



## Basic Information

**Catalog Number:**

Ag19817

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6<sup>°</sup>His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

ILEKGHAKHWFSSLRPSRNVVFSIIEHFWRELLSQPG  
AVADASVPRTRGEELTGLMAGSQKRELTARVLHQ  
GQQLDSVMRGLDKMESDLEVADRLLTELESPAWWWPF  
SSKLWKTPEETKPREDVSMTSCEPFGKEGILIKIPAVIS  
HRTESHVKPGRLTVLVSGLEIHDSSLLMHRFEREDVD  
DIKVHSPYEISIRQRFIGKPDMA YRLISAKMPEVIPLEV  
QFSKKMELLEDAVLRSARTSSPAEKSCSVWHAASGL  
MGCTLHREPPAGDQEGTALHLQTSLPALSEADTQELT  
QILRRMKGLALEAESELERQDEALDGVAAAVDRATLTI  
DKHNRMRKRLT

(115-464 aa encoded by BC018760)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.