

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

DISP2 Fusion Protein



Basic Information

Catalog Number:

Ag20127

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MDGDS S S S S G G S G P A P G P G P E G E Q R P E G E P L A P D G G
S P D S T Q T K A V A P E A S P E R S S L H S C P L E D P S S S S G P P P
T T S T L Q P V G P S S P L A P A H F T Y P R A L Q E Y Q G G S S L P G L
G D R A A L C S H G S S L S P A P S Q R D G T W K P P A V Q H H V
V S V R Q E R A F Q M P K S Y S Q L I A E W P V A V L M L C L A V I F
(1-181 aa encoded by BC151833)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.