

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

PFAS Fusion Protein



Basic Information

Catalog Number:

Ag20162

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MVRVSVNGAVVLEEPVGE LRALWEETS FQ LDR LQAE
PRCVAEEERGLRERMGPSYCLPPTFPKASVPREPGGGS
PRVA ILREEGSNGDREMADAFHLAGFEVWDVTMQDL
CSGAI GLDTRGVA FVGGFSYADVLGSAKWAAA VT
FHPRAGAE LRRFRKRPD TFS LGVCN GCQLLALLGWVG
GDPNEDAAEMGPDSQP ARPGLLRHNL SGRYESRWA
SVRVGPGPALMLRGM EGA VLPVWSAHGEGYVAFSSP
ELQAQIEARGLAPLHWADDG NPT EQYPLNPN GSPG
GVAGI CSCDGRHLA VMPHPERAVRPWQWAWRPPPF
DTLTTSPWLQLSINARNWTLEGSC
(989-1338 aa encoded by BC167158)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.