

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

PFAS Fusion Protein



Basic Information

Catalog Number:

Ag20202

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MVRVSVNGAVVLEEPVGE LRALWEETS FQLDRLQAE
PRCV AEEERGLRERMGPSYCLPPTFPKASVPREPGGPS
PRVA I LREEGSNGDREMADAFHLAGFEVWDVTMQDL
CSGAI GLDTRGVA FVGGFSYADVLGSAKWAAA VT
FHPRAGAE LRRFRKRPDTFSLGVCNGCQLLALLGWVG
GDPNEDAAEMGPDSQPAPRG LLLRHNL SGRYESRWA
SVRVGPGPALMLRGM EGA VLPVWSAHGEGYVAFSSP
ELQAQIEARGLAPLHWADDG NPT EQYPLNPN GSPG
GVAGI CSCDGRHLA VMPHPERAVRPWQWAWRPPPF
DTLTTSPWLQLSINARNW TLEGSC
(989-1338aa encoded by BC167158)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.