

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# HPCA Fusion Protein



## Basic Information

**Catalog Number:**

Ag23813

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MGKQNSKLRPEVLQDLRENTFTDHELQEWYKGF<sub>1</sub>LD  
CPTGHLTVDEFKIKIYANFFPYGDASKFAEHVFR<sub>2</sub>FTD  
NGDGTIDFREFI<sub>3</sub>I<sub>4</sub>ALS<sub>5</sub>VTS<sub>6</sub>RG<sub>7</sub>LE<sub>8</sub>Q<sub>9</sub>KL<sub>10</sub>W<sub>11</sub>A<sub>12</sub>F<sub>13</sub>S<sub>14</sub>MY<sub>15</sub>DL<sub>16</sub>G  
NGYISRSEMLEIVQAIYKMVSSVMKMPED<sub>17</sub>ST<sub>18</sub>PE<sub>19</sub>K<sub>20</sub>RT  
DKIFRQMDT<sub>21</sub>NND<sub>22</sub>G<sub>23</sub>KL<sub>24</sub>S<sub>25</sub>LEEFIRGA<sub>26</sub>KSD<sub>27</sub>PS<sub>28</sub>IV<sub>29</sub>RL<sub>30</sub>L<sub>31</sub>Q<sub>32</sub>CD<sub>33</sub>P  
SSASQF

(1-193 aa encoded by BC009846)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.