

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

LDB2 Fusion Protein



Basic Information

Catalog Number:

Ag2460

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSSTPHDPFYSSPFGPFYRRHTPYMVQPEYRIYEMNK
RLQSRTESDNLWWDAFATEFFEDDATLTLSCLEDG
PKRYTIGRTLIPRYFSTVFEGGVDLYYILKHSKESYHNS
SITVDCDQCTMVTQHKGKPMFTKVCTEGRLLLEFTFDD
LMRIKTWHFTIRQYRELVPRSILAMHAQDPQVLDQLS
KNITRMGLTNFTLNYLRCLVILEPMQELMSRHKTYNLSP
RDCLKTCLFQKWQRMVAPPAEPTRQPTTKRRKRKNS
TSSTSNSAGNNANSTGSKKTTAANLSLSSQVPDV
MVGPEPTLMGGEGFGEDEDERLITRENTQYDAANGMD
DEEDFNNSPALGNNSPWNSKPPATQETKSENPPQA
SQ

(1-373 aa encoded by BC034019)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.