

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ST7 Fusion Protein



## Basic Information

**Catalog Number:**  
Ag2547

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, T-GST, with N-terminal GST.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Validated Application:**  
Blocking peptide

**Peptide Sequence:**  
LRPLLGGVDNNSNNSDSNRQSVSECKVWRNPLNLF  
GAEYNRYTWWTGREPLTYDDMNSAQDHQTFFTCD  
DHLRPADAIMQKAWRERNPQARISAAHEALEINECA  
TAYILLAEAEATTIAEAEKLFKQALKAGDGCYRRSQ  
QLQHHGSQYEAQHRRDNTNVLVYIKRRLAMCARRL  
GRTR EAVKMMRDLMKEFPLLSMFNIHENLLELLE  
LQAYAD VQAVLAKYDDISLPKSATICYTAALLKAR  
AVSDKFSPE AASRRGLSTAEMNAVEAIHRAVEF  
NPHVPHYKYLLEMKS LILPPEHILKRGDSEAI  
YAFFHLAHWKRVEGALNLLHC TWEGTFRMIPY  
PLEKGHLYPYPICTETADRELLPSFHE VSVYP  
KKE  
(101-481 aa encoded by BC030954)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.