

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

LOXL4 Fusion Protein



Basic Information

Catalog Number:

Ag2669

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

TEAMVACRQLGLGF~~A~~I~~H~~AYKETWFW~~S~~GTPRAQE~~V~~V
MSGVRC~~S~~GTE~~L~~ALQ~~Q~~CQRHG~~P~~VHCSHG~~G~~GRFLAGV
SCMDSAPDLVMNAQLVQETAYLED~~R~~PLS~~Q~~LYCAHEE
NCLSKSADHMDWPYGYRRLLRFSTQIYNLGR~~T~~DFR~~P~~K
TGRDSWVWHQC~~H~~RHYHSIEVFTHYDLL~~T~~LNGSKVAE
GHKASFCLE~~D~~TNCPTGLQRRYACANFGEQGV~~T~~VGC
WDTYRH~~D~~ICQWVDITDVGP~~G~~NYIFQVIVNPHYEVA
ESDFSNNMLQCRC~~K~~YD~~G~~HRVWLHNCH~~T~~GNSYPANA
ELSLEQEQLRNNLI

(457 - 756 aa encoded by BC013153)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.