

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# NRP1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag2720

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

NDKCGDTIKIESPGYLTSPGYPHSYHPSEKCEWLIQAP  
DPYQRIMINFNPHFDLEDRDCKYDYVEVFDGENENG  
HFRGKFCGKIAPPPVVSQPFLLIKFVSDYETHGAGFS  
IRYEIFKRGPESQNYTTPSGVIKSPGFPEKYPNSLECT  
YIVFAPKMSEIILEFESFDLEPDSNPPGGMFCRYDRLEI  
WDGFPDVGPHIGRYCGQKTPGRIRSSSGILSMVFYTD  
SAIAKEGFSANYSVLQSSVSEDFKCMALGMESGEIH  
SDQITASSQYSTNWSAERSRLNYPENGWTPGEDSYR  
EWIQVDLGLLRFVAVGTQGAISKETKKKYYVKTYKID  
VSSNGEDWITIKENKPVLFQGNTPDVTVVAVFPK  
PLITRFVRIKPATWETGISMRFEVYGCKITDYPCSGML  
GMVSGLISDSQITSSNQGDRNWMPENIRLVTSRSGW  
ALPPAPHSYINELQIDLGEEKIVRGIIQGGKHRENKV  
FMRKFKIGYSNNGSDWKMI MDDSKRKA KSFEGNNNY  
DTPELRTFPALSTRFIRIYPERATHGGLGLRMELLGCEV  
EAPTAGPTTPNGNLVDECCDDQ  
(24-607aa encoded by BC007533)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.