

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SAV1 Fusion Protein



Basic Information

Catalog Number:

Ag3441

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MLSRKKTINEVSKPAEVQGKYVKKETSPLLRNLMPFI
RHGPTIPRRTDCLPDSSPNAFSTSGDVVSRNQSFRLR
PIQRTPHEIMRRESNRLSAPSYLARSADVPREYGSSQ
SFVTEVSFAVENGDSGSRYYSNFFDQQRKPLGDR
AHEDYRYEYNHDLFQRMPQNGRHSAGIGRVAATS
LGNLNTHGSEDLPLPGWSDWTMRGRKYIDHNTN
TTHWSHPLERGLPPGWVERVESSEFGTYVDHTNKKA
QYRHPCAPSVPRYDQPPPVTYQPQQTERNQSLLVPA
NPYHTAEIPDWLQVYARAPVKYDHIKWFELQADLD
TYQGMLKLLFMKELEQIVKMYEAYRQALLTELENRKQR
QQWYA

(1-376 aa encoded by BC020537)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.