

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# FOXI1 Fusion Protein



## Basic Information

Catalog Number:

Ag3958

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSSF~~DL~~PAPSPPRCSQ~~Q~~FPSIGQEPPEMNL~~Y~~ENFFH  
PQGV~~P~~SPQRPSFEGGGEYGATPNPVLWFNGPTMT~~PP~~  
PYLPGPNAS~~P~~FLPQAYGVQRLLPSVSLGGSDLGWL  
PIPSQEELMKLVRPPYSY~~S~~ALIAMAIHGAPDKRL~~T~~LSQI  
YQYVADNFPFY~~N~~KS~~K~~AGWQNSIRHNLSLND~~C~~FKKVP  
RDEDDPA~~V~~YVSGGSPTSHPLVTPGLSPEPSDK~~T~~GQNSL  
TFNSF~~S~~PLTNLSNHSGGDWANPMP~~T~~NMLS~~V~~GGSVL  
SQFSPHFYNSVNTSGVLYPREGTEV  
(1-283 aa encoded by BC029778)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.