

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# JAK3 Fusion Protein



## Basic Information

**Catalog Number:**  
Ag4298

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Validated Application:**  
Blocking peptide

**Peptide Sequence:**  
MAPPSEETPLIPQRSCSLLSTEAGALHVLLPARGPGPP  
QRLSFSFGDHLAEDLCVQAAKASGILPVYHSLFALAT  
EDLSCWFPPSHIFSVEDASTQVLLYRIRFYFPNWFGL  
KCHRFGLRKDLASAILDLPVLEHLFAQHRSDLVSGRLP  
VGLSLKEQGECLSLAVLDLARMAREQAQRPGELLKTV  
SYKACLPPSLRDLIQGLSFVTRRRIRRTVRRALRRVAAC  
QADRHSLMAKYIMDLERLDPAGAAETFHVGLPGALG  
GHDGLLLRVAGDGGIAWTQGEQEVLPFCDFPEIV  
DISIKQAPRVGPAGEHRLVTVTRTDNQILEAEFPGLP  
ALSFVALVDGYFR  
(1-350 aa encoded by BC028068)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.