

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# SCARF1 Fusion Protein



## Basic Information

Catalog Number:

Ag4621

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LACCACCCWAPRSDLKDRPARDGATVSRMKLQVWG  
TLTSLGSTLPCRSLSHKL PWVTVSHHDPEVPFNHSFI  
EPPSAGWATDDSFSSDPESGEADEVPA YCVPPQEG  
MVPVAQAGSSEASLAAGAFPPPEDASTPFAIPRTSSL  
ARAKRPSVSFAEGTKFAPQSRSSGELSSPLRKPRLS  
RGAQSGPEGREAEESTGPDEAEAPESFPAASP GDS  
ATGHRRPPLGSRVVAEHVEAIEG SVQESSGPVTIYM  
LAGKPRGSEGPVRSVFRHFGSFQKGQAEAKVKRAIPK  
PPRQALNRKKGSPGLASG SVGQSPNSAPKAGLPGAT  
GPM AVRPEEAVRGLGAGTESSRRAQEPVSGCGSPEQ  
DPQKQAE EERQEEPEYENVVPI SRPPEP  
(438-830 aa encoded by BC039735)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.