

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

KERA Fusion Protein



Basic Information

Catalog Number:

Ag4627

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

RQVYEVHSDSDWTIHD F E C P M E C F C P P S F P T A L Y C E
N R G L K E I P A I P S R I W Y L Y L Q N N L I E T I P E K P F E N A T Q L R
W I N L N K N K I T N Y G I E K G A L S Q L K K L L F L F L E D N E L E E V P
S P L P R S L E Q L Q L A R N K V S R I P Q G T F S N L E N L T L L D L Q N
N K L V D N A F Q R D T F K G L K N L M Q L N M A K N A L R N M P P R L
P A N T M Q L F L D N N S I E G I P E N Y F N V I P K V A F L R L N H N K L
S D E G L P S R G F D V S S I L D L Q L S H N Q L T K V P R I S A H L Q H L
H L D H N K I K S V N V S V I C P S P S M L P A E R D S F S Y G P H L R Y L
R L D G N E

(24-331aa encoded by BC032667)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.