

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SENP2 Fusion Protein



Basic Information

Catalog Number:

Ag4664

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MKETIQGTGSWGPEPPGPIPPAYSSPRRERLRWPPP
PKPRLKSGGGFGPDGSGTTVPARRLPVPRPSFDASA
SEEEEEEEEEDEEEVAAWRLPPRWSQLGTSQRP
RPSRPTHKTC SQRRRRAMRAFRMLLYSKSTSLTFHW
KLWGRHRGRRRGLAHPKNHLSPPQGGATPQVPSPCC
RFDSPRGPPPPRLGALLGALMAEDGVRGSPVPVSGPPM
EEDGLRWTPKSPDPDGLSCTLPNGFGGQSGPEGE
RSLAPPDASILSNVCSIGDHVAQELFQGSGLGMAEE
AERPGEKAGQHSPREEHVTCVQSILDEFLQTYGSLIP
LSTDEVVEKLEDIFQQEFSTPSRKLGLVLQLIQSYQRMP
GNAMVRGFRVAYKRHLVLTMDLGLTYGQNLNDQV
MNMVYDGLVMDTVPEKVHFFNSFFYDKLRTKGYDGVK
RWTKNVDIFNKELLIPHLEVHWSLISVDVRRRTITYF
DSQRTLNRRCPIKIAKYLQAEAVKKDRDLFHQGWKG
YFKMNVARQNNDSDCGAFVLQYCKHLALSQPF5FT
QQDMPKLRRIYKELCHCKLTV
(1-589aa encoded by BC080658)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.