

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# FCRL1 Fusion Protein



## Basic Information

Catalog Number:

Ag4724

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6<sup>+</sup>His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MLPRLLLIICAPLCEPAELFIASPSHPTEGSPVTLTCKM  
PFLQSSDAQFQFCFFRDTRALGPGWSSSPKLQIAAM  
WKEDTGSYWCEAQTMASKVLRSSQINVHRVPVA  
DVSLETQPPGGQVMEGDRLVLI CSVAMGTGDITFLW  
YKGA VGLNLQSIKTQRSLTAEYEI PSVRESDAEQYYCV  
AENGYGPSGLVSI TVRIPVSRPILMLRAPRAQA AVE  
DVLELHCEALRGSPPIYWFYHEDITLGSRSAPSGGGA  
SFNLSLTEEHSNGYSCEANGLGAQRSEAVTLNFTVP  
TGARSNHLTSGVIEGLLSTLGPATVALLFCYGLKRKIGR  
RSARDPLRSLPPLQEFYTLNSPTPGQLQPIYENVNV  
VSGDEVYSLAYNQPEQESVAAETLGHMEDKVSLEI  
YSRLRKANITDVDYEDAM

(1-429aa encoded by BC033690)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.