

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

DARS2 Fusion Protein



Basic Information

Catalog Number:

Ag4809

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LI E G L L Q Y S W P N D K D P V V V P F P T M T F A E V L A T Y G T D K
P D T R F G M K I I D I S D V F R N T E I G F L Q D A L S K P H G T V K A I
C I P E G A K Y L K R K D I E S I R N F A A D H F N Q E I L P V F L N A N R
N W N S P V A N F I M E S Q R L E L I R L M E T Q E E D V V L L T A G E H
N K A C S L L G K L R L E C A D L L E T R G V V L R D P T L F S F L W V V D
F P L F L P K E E N P R E L E S A H H P F T A P H P S D I H L L Y T E P K K A
R S Q H Y D L V L N G N E I G G G S I R I H N A E L Q R Y I L A T L L K E D V
K M L S H L L Q A L D Y G A P P H G G I A L G L D R L I C L V T G S P S I R D
V I A F P K S F R G H D L M S N T P D S V P P E E L K P Y H I R V S K P T D
S K A E R A H

(296-645 aa encoded by BC045173)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.