

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

PLA1A Fusion Protein



Basic Information

Catalog Number:

Ag5467

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MPPGPWESC FVWGG L L W L S V G S S G D A P P T P Q P K C A
D F Q S A N L F E G T D L K V Q F L L F V P S N P S C G Q L V E G S S D L
Q N S G F N A T L G T K L I I H G F R V L G T K P S W I D T F I R T L L R A T
N A N V I A V D W I Y G S T G V Y F S A V K N V L G V S E S S I H I I G V
S L G A H V G G M V G Q L F G G Q L G Q I T G L D P A G P E Y T R A S V
E E R L D A G D A L F V E A I H T D T D N L G I R I P V G H V D Y F V N G
G Q D Q P G C P T F F Y A G Y S Y L I C D H M R A V H L Y S A L E N S C
P L M A F P C A S Y K A F L A G R C L D C F N P F L L S C P R I G L V E Q G
G V K I E P L P K E V K V Y L L T S S A P Y C M H H S L V E F H L K E L R N
K D T N I E V T F L S S N I T S S S K I T I P K Q Q R Y G K G I A H A T P Q
C Q I N Q V K F K F Q S S N R V W K K D R T T I I G K F C T A L L P V N D
R E K M V C L P E P V N L Q A S V T V S C D L K I A C V
(1-440aa encoded by BC047703)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.