

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# KRT36 Fusion Protein



## Basic Information

**Catalog Number:**

Ag5483

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, T-HIS, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MLRSGLSGLGSLPGSYLSSECHTSGFVGGWFCEG  
SFNGSEKETMQFLNDRLANYLEKVRQLERENAELESRI  
QEWYEFQIPYICPDYQSYFKTIEDFQQKILLTKSENAR  
LVLQIDNAKLAADDFTKYETELSLRQLVEADINGLRR  
LDELTLCKADLEAQVESLKEELMCLKKNHEEEVSVLRC  
QLGDRLNVEVDAAPPVLDNKILEDMRCQYEALVENN  
RRDVEAWFNTQTEELNQVVSSSEQLQCCQTEIIEI  
RRTVNALEIELQAQHSMRNSLESTLAETEARYSSQLA  
QMQCLISNVEAQLSEIRCDLERQNQEYQVLLDVKARL  
EGEIATYRHLEGEDCKLPPQPCATACKPVIRVPSVPP  
VPCVPSVPCTPAPQVGTQIRITITEIRDGKVISSREHV  
QSRPL

(1-417 aa encoded by BC043581)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.