

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

POLK Fusion Protein



Basic Information

Catalog Number:

Ag6194

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MDSTKEKCD^{SY}KD^{LL}RLMGLNDNKA^GMEGLDKEKIN
KIIMEATK^GSRFYGNELKKEKQV^NQRIENMMQ^KKAQI
TSQQLRKAQLQVDRFAMELEQ^SSRNLSNTIVHIDMDA
FYAAVEMRDNP^ELKDKPIAVGSM^SMLSTSNYHARRFG
VRAAMP^GFI^AKRLCPQLIIVPPNFDKYRAV^SKEVKEIL
ADYDPNFMAMSLDEAYLNITKHLEERQ^NWPEDKRRYF
IKMGSSVENDNPGKEV^NKLSEHERSISPL^LFEESPSDV
QPPGDPFQVNFEEQ^NNPQILQNSVVF^GTSAQEVVK
EIRFRIEQK^TTLTASAGIAPNTMLAKV^CSDKNKPNGQY
QILPNRQAV^MDFIK

(1-347 aa encoded by BC050718)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.