

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PDHB Fusion Protein



## Basic Information

**Catalog Number:**

Ag6857

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MAAVSGLVRRPLREVSGLLKRRFHWTAPAALQVTVR  
DAINQGMDEELERDEKVFLLGEEVAQYDGA YKVS RGL  
WKYGDKR I IDTPISEMFGAGI AVGAAMAGLRPI CEF  
MTFNFSMQAIDQVINSAAKTYMSSGGLQPVPIVFRG  
PNGASAGVAAQHSQCFAAWYGHCPGLKVVSPWNS  
EDAKGLIKSAIRDNNP VVLENELMYGV PFEFPPEAQ  
SKDFLIPIGKAKIERQGTHTV VSHSRPVGHCL EAAAV  
LSKEGVECEVINMRTIRPMDMETIEASVMKTNHLVTV  
EGGWPQFVGAEICARIMEGPAFNFLDAPAVRVTGA  
DVPMPYAKILEDNSIPQVKDIIFAIKKTLNI  
(1-359 aa encoded by BC000439)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.