

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# HNRNPF Fusion Protein



## Basic Information

**Catalog Number:**

Ag7128

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MMLGPEGGEGFVVKLRGLPWSCSVEDVQNFLSDCTI  
HDGAAGVHFIYTRGRQSGEAFVELGSEDDVKMALK  
KDRESMGHRYIEVFKSHRTEMWVWKHSGPNSADSA  
NDGFVRLRGLPFGCTKEEIVQFFSGLEIVPNGITLPVD  
PEGKITGEAFVQFASQELAEKALGKHKERIGHRYIEVF  
KSSQEEVRSYSDPPLKFMSVQRPGPYDRPGTARRYIGI  
VKQAGLERMRPGAYSTGYGGYEEYSGLSGDYGFDDL  
FGRDLSYCLSGMYDHRVGDSEFTVQSTTGHCVMHMRG  
LPYKATENDIYNFFSPLNPRVHIEIGPDGRVTGEADV  
EFATHEEAVAAMSKDRANMQHRYIELFLNSTTGASN  
GAYSSQVMQGMGVSAAQATYSGLESQSVSGCYGA  
GYSQQNSMGGYD

(1-415 aa encoded by BC001432)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.