

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ATG4B Fusion Protein



## Basic Information

Catalog Number:

Ag7261

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MDAATLTYDTRLRFAEFEDFPETSEPVWILGRKYSIFTE  
KDEILSDVASRLWFTYRKNFPAIGGTGPTSDTGWGCM  
LRCGQMIFAQALVCRHLGRDWRWTQRKRQPDYSYFS  
VLNAFIDRKDSYYSIHQIAQMGVGEKSIQWYGP  
TVAQVLKLVFDTWSSSLAVHIAMDNTVVMEEIRRL  
CRTSVPCAGATAFPADSDRHCNGFPAGAEVTNRPS  
WRPLVLLIPLRLGLTDINEAYVETLKHCFMMPQSLGVI  
GGKPNSAHYFIGYVGEELIYDPHTTQPAVEPTDGCFI  
PDEFHQCQHPPCRMSIAELDPSIAGVGFCKTEDDFND  
WCQQVKKLSLLGGALPMFELVEQQPSHLACPDVNL  
SLDSSDVERLERFFDSEDEDFEILSL  
(1-393 aa encoded by BC000719)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.