

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

ANT3 Fusion Protein



Basic Information

Catalog Number:

Ag7703

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MTEQAI^SFAKDFLAGGIAAAISKTAVAPIERVKLLQV
QHASKQIAADKQYKGI^VDCIVRIPKEQGVLSFWRGNL
ANVIRYFPTQALNFAFKDKYKQIFLGGVDKHTQFWRY
FAGNLASGGAAGATSLCFVYPLDFARTRLAADV^GKS
GTEREFRGLGDCLVKITKSDGIRGLYQGFVSVQGIIY
RAAYFGVYDTAKGMLPDPKINTHIVVSWMIAQTVTAV
AGVVSYPFDTVRRRMMMQSGRKGADIMYGTVD^CW
RKIFRDEGGKAFFKGAWSNVL^RGMGGAFVLVLYDELK
KVI

(1-298aa encoded by BC008935)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.