

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SPTLC2 Fusion Protein



Basic Information

Catalog Number:

Ag7774

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

CLILSDELNHASLVLGARLSGATIRIFKHNNMQSLEKLL
KDAIVYQGQPRTRRPWKKILILVEGIYSMEGSIVRLPEVI
ALKKKYKAYLYLDEAHSIGALGPTGRGVVEYFGLDPED
VDVMMGTFTKSF GASGGYIGGKELIDYLRTHSHSAV
YATSLSPVVEQIITSMKCIMGQDGTSLGKECVQQLA
ENTRYFRRLKEMGFIIYGNEDSPVVPLMLYMPAKIGA
FGREMLKRNI GVVVVGFPATPIIESRARFCLSAHTKE
ILD TALK EID E V G D L L Q L K Y S R H R L V P L L D R P F D E T T Y E E
T E D

(254-562 aa encoded by BC005123)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.