

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ABI3 Fusion Protein



## Basic Information

**Catalog Number:**

Ag8164

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MAELQQLQEF E I P T G R E A L R G N H S A L L R V A D Y C E D N Y  
V Q A T D K R K A L E E T M A F T T Q A L A S V A Y Q V G N L A G H T L  
R M L D L Q G A A L R Q V E A R V S T L G Q M V N M H M E K V A R R E  
I G T L A T V Q R L P P G Q K V I A P E N L P L T P Y C R R P L N F G C L  
D D I G H G I K D L S T Q L S R T G T L S R K S I K A P A T P A S A T L G R P  
P R I P E P V H L P V V P D G R L S A A S S A S L A S A G S A E G V G G  
A P T P K G Q A A P P A P L P S S L D P P P P P A A V E V F Q R P P T L  
E E L S P P P D E E L P L D L P P P P L D G D E L G L P P P P G F G P  
D E P S W W P A S Y L E K V V T L Y P Y T S Q K D N E L S F S E G T V I C  
V T R R Y S D G W C E G V S S E G T G F F P G N Y V E P S C  
(1-366 aa encoded by BC007780)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.