

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

MARCO Fusion Protein



Basic Information

Catalog Number:

Ag8219

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

PGPPAEKGA KGAMGRDGATG P SGPQ GPPGVKGEAG
LQGPQ G A P G K Q G A T G T P G P Q G E K G S K G D G G L I G P K G
E T G T K G E K D L G L P G S K G D R G M K G D A G V M G P P G A Q
G S K G D F G R P G P P L A G F P G A K G D Q Q P G L Q G V P G P
P G A V G H P G A K G E P G S A G S P G R A G L P G S P G S P G A T G L
K G S K G D T G L Q G Q G R K G E S G V P G P A G V K G E Q G S P G
L A G P K G A P G Q A G Q K G D Q G V K G S S G E Q G V K G E K G E R
G E N S V S V R I V G S S N R G R A E V Y Y S G T W G T I C D D E W Q N
S D A I V F C R M L G Y S K G R A L Y K V G A G T G Q I W L D N V Q C R
G T E S T L W S C T K N S W G H H D C S H E E D A G V E C S V
(171-520 aa encoded by BC016004)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.