

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# AGXT2L2 Fusion Protein



## Basic Information

Catalog Number:

Ag8374

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MGKSI~~NG~~HPVACVAATQPVARAFEATGVEYFNTF  
GGSPVSCAVGLAVLN~~VLE~~KEQLQDHATSVG~~S~~FLMQL  
LGQQKIKHPIVGDVRGVGLFIV~~DL~~IKDEATRTPATEE  
AA~~Y~~LVSR~~L~~KENYVLLSTDG~~P~~GRN~~L~~KFKPPMCFSLDNA  
RQV~~V~~AKLDA~~I~~LTD~~M~~E~~E~~KV~~R~~SC~~E~~T~~L~~R~~L~~Q~~P~~  
(276-450 aa encoded by BC008009)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.