

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SSR3 Fusion Protein



Basic Information

Catalog Number:

Ag8699

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LVAFA YK N V K F V L K H K V A Q K R E D A V S K E V T R K L S E A D
N R K M S R K E K D E R I L W K K N E V A D Y E A T T F S I F Y N N T L F L
V V V I V A S F F I L K N F N P T V N Y I L S I S A S S G L I A L L S T G S K
(72-185aa encoded by BC017203)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.