

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# SERPINB2 Fusion Protein



## Basic Information

**Catalog Number:**

Ag9319

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MEDLCVANTLFA~~N~~LFKHLAKASPTQNLFLSPW~~S~~ISST  
MAMVYMGSRG~~S~~TE~~D~~QMAKVLQFNEVGANAVTPMT  
PENFTSCGF~~M~~Q~~Q~~I~~Q~~KGSYPDA~~I~~LQAQAADKIHSSFR  
SLSSAINASTGNYLLESV~~N~~KLFGEKSASF~~R~~E~~E~~YIRLCQK  
YYSSE~~P~~QAVDFLECAEEARKKIYSWVITQTKGKIPNLL  
PEGSVDGDTRMVLVNAVYFKGKWKTPFEK~~K~~LNGLYPF  
RVNSAQRTPVQMMYLREKLNIGYIEDLKAQILELPYAG  
DVSMFLLLPDEIADVSTGLELLESEITYDKLNKWT~~S~~KDK  
MAEDEV~~E~~VYI~~P~~QFKLEE~~H~~YELRSILRSMGMEDAFNKG  
RANFSGMSERNDLFLSEV~~F~~HQAMVDVNEEGTEAAA  
GTGGVMTGR~~T~~G~~H~~GGPQFVADHPFL~~L~~MHKITNCILF  
FGRFSSP

(1-415 aa encoded by BC012609)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.