

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ELP3 Fusion Protein



## Basic Information

Catalog Number:

Ag9993

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MRQKRKGDLSPAELMMLTIGDVIKQLIEAHEQGKDIDL  
NKVKTATAAKYGLSAQPRLVDIIAAVPPQYRKVLMPKL  
KAKPIRTASGIAVVAVMCKPHRCPHISFTGNICVYCP  
GGPDSDFEYSTQSYTGYEPTSMRAIRARYDPFLQTRH  
RIEQLKQLGHSVDKVEFIVMGGTFMALPEEYRDYFIRN  
LHDALSGHTSNNIYEAVKYSERSLTKCIGITIEIRPDYC  
MKRHLSMMLTYGCTREIGVQSVYEDVARDTNRGHT  
VKAVCESFHAKDSGFKVVAHMMPDLPNVGLERDIE  
QFTEFFENPAFRPDGLKLYPTLVIRGTGLYELWKSGRY  
KSYSPSDLVELVARILALV

(1-356 aa encoded by BC001240)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.