

For Research Use Only

CoraLite® Plus 488-conjugated LSM12 Monoclonal antibody



Catalog Number:CL488-67029

Basic Information

| | | |
|---|--|---|
| Catalog Number: CL488-67029 | GenBank Accession Number: BC044587 | Purification Method: Protein A purification |
| Size: 100ul , Concentration: 1000 µg/ml by Nanodrop; | GeneID (NCBI): 124801 | CloneNo.: 2E3G3 |
| Source: Mouse | Full Name: LSM12 homolog (S. cerevisiae) | Excitation/Emission maxima wavelengths: 488 nm / 515 nm |
| Isotype: IgG2a | Observed MW: 22 kDa | |
| Immunogen Catalog Number: AG28782 | | |

Applications

Tested Applications:
FC (Intra)

Species Specificity:
Human, Mouse, Rat

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human LSM12 (CL488-67029, Clone:2E3G3) (red), or 0.4 ug Mouse IgG2a Isotype Control (CL488-66360-2, Clone: K11A1B2A2) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).