

For Research Use Only

CoraLite® Plus 488-conjugated MESDC2 Monoclonal antibody



Catalog Number: **CL488-67441**

Basic Information

Catalog Number: CL488-67441	GenBank Accession Number: BC009210	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 23184	CloneNo.: 1A2E7
Source: Mouse	Full Name: mesoderm development candidate 2	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG2a	Calculated MW: 26 kDa	Excitation/Emission maxima wavelengths: 488 nm / 515 nm
Immunogen Catalog Number: AG1403		

Applications

Tested Applications: IF	Positive Controls: IF : HEK-293 cells,
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using CoraLite® Plus 488 MESDC2 antibody (CL488-67441, Clone: 1A2E7) at dilution of 1:100.