

For Research Use Only

CoraLite® Plus 488-conjugated PDS5A Monoclonal antibody



Catalog Number:CL488-67520

Basic Information

Catalog Number: CL488-67520	GenBank Accession Number: BC009650	Purification Method: Protein G purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 23244	CloneNo.: 1C8D4
Source: Mouse	Full Name: PDS5, regulator of cohesion maintenance, homolog A (S. cerevisiae)	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG1	Calculated MW: 1297 aa, 147 kDa	Excitation/Emission maxima wavelengths: 488 nm / 515 nm
Immunogen Catalog Number: AG11525		

Applications

Tested Applications: IF	Positive Controls: IF : SMMC-7721 cells,
Species Specificity: Human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed SMMC-7721 cells using CoraLite® Plus 488 PDS5A antibody (CL488-67520, Clone: 1C8D4) at dilution of 1:200.