

For Research Use Only

CoraLite® Plus 488-conjugated FARSB Monoclonal antibody



Catalog Number: CL488-67924

Basic Information

Catalog Number: CL488-67924	GenBank Accession Number: BC017783	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 10056	CloneNo.: 1H10C2
Source: Mouse	Full Name: phenylalanyl-tRNA synthetase, beta subunit	Recommended Dilutions: IF 1:50-1:500
Isotype: IgG2a	Calculated MW: 589 aa, 66 kDa	
Immunogen Catalog Number: AG9695	Observed MW: 66 kDa	

Applications

Tested Applications: IF	Positive Controls: IF : U-251 cells,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed U-251 cells using CoraLite® Plus 488 FARSB antibody (CL488-67924, Clone: 1H10C2) at dilution of 1:200.



Immunofluorescent analysis of (-20°C Methanol) fixed U-251 cells using CoraLite® Plus 488 FARSB antibody (CL488-67924, Clone: 1H10C2) at dilution of 1:200.