

For Research Use Only

CoraLite® Plus 750 Anti-Human CD1a (HI149)



Catalog Number: CL750-65172

Basic Information

Catalog Number: CL750-65172	GenBank Accession Number: BC031645	Purification Method: Affinity purification
Size: 100tests , 5 µl/test	GeneID (NCBI): 909	CloneNo.: HI149
Source: Mouse	Full Name: CD1a molecule	Excitation/Emission maxima wavelengths: 755 nm / 780 nm
Isotype: IgG1, kappa	Calculated MW: 327 aa, 37 kDa	

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ MOLT-4 cells were surface stained with 5 ul CoraLite® Plus 750 Anti-Human CD1a (CL750-65172, Clone:HI149) or Mouse IgG1 Isotype Control. Cells were not fixed.



1X10⁶ MOLT-4 cells were surface stained with 5 ul CoraLite® Plus 750 Anti-Human CD1a (CL750-65172, Clone:HI149) or Mouse IgG1 Isotype Control. Cells were not fixed.