

For Research Use Only

FITC Plus Anti-Human CD58 (TS2/9)

Catalog Number: FITC-65193



Basic Information

Catalog Number: FITC-65193	GenBank Accession Number: BC005930	Purification Method: Protein G purification
Size: 100tests , 5 µl/test	GeneID (NCBI): 965	CloneNo.: TS2/9
Source: Mouse	Full Name: CD58 molecule	Excitation/Emission maxima wavelengths: 495 nm / 524 nm
Isotype: IgG1, kappa	Calculated MW: 28 kDa	

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface stained with 5 ul FITC Anti-Human CD58 (FITC-65193, Clone:TS2/9) or unstained. Cells were not fixed. Monocytes were gated.