

For Research Use Only

Calpain 10 Polyclonal antibody

Catalog Number: 10419-1-AP



Basic Information

Catalog Number:

10419-1-AP

Size:

150ul , Concentration: 350 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0661

GenBank Accession Number:

BC004260

GeneID (NCBI):

11132

UNIPROT ID:

Q9HC96

Full Name:

calpain 10

Calculated MW:

75 kDa

Observed MW:

57 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Jurkat cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 10419-1-AP (Calpain 10 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.