

For Research Use Only

L3MBTL3 Polyclonal antibody

Catalog Number: 14070-1-AP



Basic Information

Catalog Number:

14070-1-AP

Size:

150ul , Concentration: 500 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5130

GenBank Accession Number:

BC060845

GeneID (NCBI):

84456

UNIPROT ID:

Q96JM7

Full Name:

l(3)mbt-like 3 (Drosophila)

Calculated MW:

88 aa, 3 kDa

Observed MW:

90 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : Jurkat cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 14070-1-AP (L3MBTL3 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.