

For Research Use Only

TRMT10A Polyclonal antibody

Catalog Number: 17294-1-AP



Basic Information

Catalog Number: 17294-1-AP	GenBank Accession Number: BC028373	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 700 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 93587	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	Full Name: RNA (guanine-9-) methyltransferase domain containing 2	
Isotype: IgG	Calculated MW: 339 aa, 40 kDa	
Immunogen Catalog Number: AG11172	Observed MW: 40 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HepG2 cells, Jurkat cells, NIH/3T3 cells
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17294-1-AP (TRMT10A antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.