

For Research Use Only

NUDT15 Polyclonal antibody

Catalog Number: 25914-1-AP



Basic Information

Catalog Number: 25914-1-AP	GenBank Accession Number: BC133015	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 µg/ml by Nanodrop;	GeneID (NCBI): 55270	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9NV35	
Isotype: IgG	Full Name: nudix (nucleoside diphosphate linked moiety X)-type motif 15	
Immunogen Catalog Number: AG23242	Observed MW: 19 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells, HeLa cells
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat lysates were subjected to SDS PAGE followed by western blot with 25914-1-AP (NUDT15 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.