

For Research Use Only

ZBED1 Polyclonal antibody

Catalog Number: 26140-1-AP



Basic Information

Catalog Number: 26140-1-AP	GenBank Accession Number: BC015030	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 400 ug/ml by Nanodrop and 267 ug/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 9189	Recommended Dilutions: WB 1:1000-1:6000
Source: Rabbit	UNIPROT ID: O96006	
Isotype: IgG	Full Name: zinc finger, BED-type containing 1	
Immunogen Catalog Number: AG23578	Observed MW: 78 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, mouse heart tissue, Jurkat cells
Species Specificity: human, mouse	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26140-1-AP (ZBED1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.