

For Research Use Only

DUB3 Polyclonal antibody

Catalog Number: 26143-1-AP

Featured Product

5 Publications



Basic Information

Catalog Number:

26143-1-AP

Size:

150ul, Concentration: 280 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23567

GenBank Accession Number:

BC157073

GeneID (NCBI):

377630

UNIPROT ID:

Q6R6M4

Full Name:

deubiquitinating enzyme 3

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : SGC-7901 cells,

Notable Publications

Author	Pubmed ID	Journal	Application
Haile Zhao	34658851	Front Pharmacol	WB,IP
Zhitao Jiang	35849971	Phytomedicine	WB
Yang Su	36040642	Breast Cancer Res Treat	WB,IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26143-1-AP (DUB3 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.