

For Research Use Only

PEX26 Polyclonal antibody

Catalog Number: 27472-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

27472-1-AP

Size:

150ul, Concentration: 600 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26529

GenBank Accession Number:

BC016280

GeneID (NCBI):

55670

UNIPROT ID:

Q7Z412

Full Name:

peroxisomal biogenesis factor 26

Calculated MW:

34 kDa

Observed MW:

30-34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : A431 cells, HEK-293 cells

Notable Publications

Author	Pubmed ID	Journal	Application
Richard G Lee	38129691	Nat Cell Biol	WB
Bianbian Yan	37957408	Dig Dis Sci	IHC, WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27472-1-AP (PEX26 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.