

For Research Use Only

ARHGAP18 Polyclonal antibody

Catalog Number: 28409-1-AP



Basic Information

Catalog Number: 28409-1-AP	GenBank Accession Number: BC107416	Purification Method: Antigen affinity purification
Size: 150ul, Concentration: 800 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 93663	Recommended Dilutions: WB 1:500-1:3000 IF/ICC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q8N392	
Isotype: IgG	Full Name: Rho GTPase activating protein 18	
Immunogen Catalog Number: AG29144	Calculated MW: 663 aa, 75 kDa	
	Observed MW: 75 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls: WB : HeLa cells, MDA-MB-453s cells IF/ICC : HeLa cells,
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28409-1-AP (ARHGAP18 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using 28409-1-AP (ARHGAP18 antibody) at dilution of 1:50 and CoraLite488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).